

	Action	Description
Registration	AiCure Registration-Patient Connect App   	<p>Register (aka activate) the participant in the AiCure Site Dashboard in-clinic at the start of the Screening or Treatment Period.</p> <ul style="list-style-type: none"> • Select participant device: an AiCure Device or the participant's personal device (BYOD). <ul style="list-style-type: none"> • If using a participant's personal device, download the Patient Connect by AiCure app in the Apple App Store or Google Play with the purple icon. • If using an AiCure Device, send a charger and case home with the participant. • Power on device and connect to WiFi. • Remind participants to keep the device powered on and charged throughout the study. We recommend connecting to WiFi whenever possible to help with data transfer. • Register the participant in AiCure.
Screening and Treatment Periods	Monitoring	<p>Monitor participants progress in the AiCure Site Dashboard throughout the study. Site coordinators should review participant medication adherence, treatment period start & end dates, and participants visit dates.</p> <ul style="list-style-type: none"> • Monitor medication adherence in the AiCure Site Dashboard. If participant non-adherence is observed (i.e. missed doses or adherence <80%), please follow up with the participant. <ul style="list-style-type: none"> ◦ Goal: Support the participant in maintaining a compliant dosing routine and inquire if there is anything preventing them from dosing. Re-educate participants on proper dosing steps with AiCure as needed. • Reclassify missed doses as 'site-confirmed' or 'not expected' in the AiCure Site Dashboard as applicable. • As applicable, please modify the participants treatment regimen in the AiCure Site Dashboard to ensure they are receiving the correct dosing instructions in the app. • The participants treatment end date is auto populated at registration to match the default protocol duration. Site coordinators can extend treatment end dates as needed in the AiCure Site Dashboard. • Site visits will automatically populate per protocol at registration. Site coordinators should modify visit dates in the AiCure Site Dashboard as needed. NOTE: Visits must be updated at least 48 hours in advanced for the in-app reminders to trigger correctly.
End of Treatment	End Participation	<p>Site coordinators must manually deactivate participants who have completed the study, early terminated, or screen failed. NOTE: The participants device must be connected to the Internet and the app opened to fully deactivate.</p> <p>Device Reset</p> <ul style="list-style-type: none"> • Confirm the participants status is 'inactive' in the AiCure Site Dashboard. • Personal Device - confirm there is no pending data for the participant on the site coordinator dashboard, then delete app from phone. • AiCure Device - Collect AiCure provisioned device and charger. Reset AiCure provisioned device for next patient.

Need help or have questions?

Visit the [AiCure Help Center](#) to submit any type of request, view articles, live chat, and more!