

How to Change the Treatment Regimen on the Dashboard

Q. How do I change the treatment regimen for my participant?

A: If your trial is configured with different treatment regimens, you can follow the steps below to change it:

Part 1: On the Dashboard

1. Log into the AiCure portal (login.aicure.com) and select your study
2. Click on the Participant ID
3. Click on the **PARTICIPANT SETTINGS** tab
4. Click on the **MANAGE PARTICIPANT** dropdown list
5. Click on **CHANGE TREATMENT REGIMEN**
6. Select the desired treatment protocol
7. Set the dosing Reminder Time
8. Select a date when the new regimen should be applied
9. Review your changes. If all is correct, click **SAVE**

An update will then be sent to the participant's device

Please Note

If the participant continues to use the app to dose without being connected to the internet to download this update, the doses may be recorded as being taken with the current regimen, and not the new one. So it is important that they update their app on the date you specified in Step 8.

PARTICIPANTS

SHOW: all completed/removed not activated
 dose in clinic today needs follow up
 active ending soon

Add Participant

Medication Adherence (AiCure confirms medication intake)	Total Number Registered	
92%	1 active	2 completed

PARTICIPANT #	START DATE	7 DAY HISTORY							TODAY (WED)	ADHERENCE
		WED	THU	FRI	SAT	SUN	MON	TUE		
001-TEST	2024-12-18	■	■	■	■	■	■	■	■	90%

SUMMARY | CALENDAR | PARTICIPANT SETTINGS

CONTACT PARTICIPANT

Treatment & Device Settings

MANAGE PARTICIPANT

CHANGE PERIOD END DATE

CHANGE TREATMENT REGIMEN

CHANGE PARTICIPANT LANGUAGE

CHANGE PARTICIPANT PHONE NUMBER

CHANGE DOSE DATE

START NEXT PERIOD

END PARTICIPATION

NEW DEVICE

CHANGE TREATMENT REGIMEN

Current Treatment Regimen
Full Step Treatment (5 Step Dosing Checks)

Select Treatment Regimen

Full Step Treatment (5 Step Dosing Checks)

Reduced Step Dosing Treatment (3 Step Dosing Treatment)

Sublingual (Sublingual)

Inhaler - Daily Dose (Inhaler modality for every day dosing)

Blister Pack & Food Reminder (Blister Pack and Food Reminders)

Care Giver Daily Dosing (Daily dosing reminder for caregivers)

Abbvie screening (Abbvie screening)

NEXT >

CHANGE TREATMENT REGIMEN

Set a reminder time for the medication

Full Step Treatment

09:00

Participants can change their reminder time from their app later

< BACK NEXT >

CHANGE TREATMENT REGIMEN

Schedule Treatment Update

Choose a date to apply treatment update:

January 2025

Su	Mo	Tu	We	Th	Fr	Sa
			1	2	3	4
5	6	7	8	9	10	11
12	13	14	15	16	17	18
19	20	21	22	23	24	25
26	27	28	29	30	31	

Date applied:

2025-01-09

< BACK NEXT >

CHANGE TREATMENT REGIMEN

Treatment Summary

Period	AiCure Treatment Period
Language	English (U.S.)
Treatment Regimen	Full Step Treatment
Reminder Times	09:00
Period Start Date	2025-01-06
Period End Date	2025-01-29
Period Duration	24 days
Apply Treatment Update	2025-01-09

< BACK SAVE

CHANGE TREATMENT REGIMEN

UPDATE PENDING

102-1234


Treatment regimen updated to Full Step Treatment for participant

Updates will take effect when the participant opens the AiCure app. No further edits can be made until update is applied to the mobile device.

OK

Part 2: On the Participant's Device

1. Instruct the participant to open the app while connected to the internet to download the update (if the app is already open, instruct them to close and reopen it)
2. Please wait while the update downloads
3. Once the update completes, the app will confirm there was an update made and the participant can proceed to using the app as normal



- Q. I am attempting to change the regimen, but the one I'm looking for is not listed. What can I do?**
- A.** The regimens listed in Part 1 Step 6 are preconfigured based on what is approved by the trial sponsors. If you cannot find the regimen you are looking for, we recommend reaching out to your sponsor team for further guidance.